

NDUFA6 Rabbit Polyclonal Antibody



CAB9831

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

13kDa

Calculated MW:

17kDa

Applications:

WB IHC

Reactivity:

Human, Mouse, Rat

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000 IHC 1:50
- 1:100

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

4700

Uniprot

P56556

Synonyms:

NDUFA6; B14; CI-B14; LYRM6; NADHB14


Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 35-154 of human NDUFA6 (NP_002481.2).


Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.


Product Images



Western blot analysis of extracts of various cell lines, using NDUFA6 antibody (CAB9831) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (CABS014) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit (CABM00020). Exposure time: 90s.



Immunohistochemistry of paraffin-embedded mouse lung using NDUFA6 antibody (CAB9831) at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded human stomach using NDUFA6 antibody (CAB9831) at dilution of 1:100 (40x lens).