

Product Information

Product SKU:	CAB9866	Gene ID:	11177	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	250kDa	Conjugate:	Unconjugated
Calculated MW:	179kDa	Isotype:	IgG

Immunogen Information

Background: The BAZ1A gene encodes the accessory subunit of the ATP-dependent chromatin assembly factor (ACF), a member of the ISWI ('imitation switch') family of chromatin remodeling complexes (summarized by Racki et al., 2009 [PubMed 20033039]).

Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	ACF1; WALp1; hACF1; WCRF180; BAZ1A
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 1347-1556 of human BAZ1A (NP_038476.2).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.