

DSTYK Rabbit Polyclonal Antibody



CAB9870

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

105kDa

Calculated MW:

20kDa/66kDa/99kDa/105kDa

Applications:

WB

Reactivity:

Human, Mouse

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

This gene encodes a dual serine/threonine and tyrosine protein kinase which is expressed in multiple tissues. It is thought to function as a regulator of cell death. Multiple transcript variants encoding different isoforms have been found for this gene.

Immunogen information

Gene ID:

25778

Uniprot

Q6XUX3

Synonyms:

DSTYK; CAKUT1; DustyPK; HDCMD38P; RIP5; RIPK5; SPG23

Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 635-884 of human DSTYK (NP_955749.1).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
