

Product Information

Product SKU:	CAB9876	Gene ID:	54538	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human

Additional Information

Observed MW:	130kDa	Conjugate:	Unconjugated
Calculated MW:	107kDa	Isotype:	IgG

Immunogen Information

Background: Predicted to enable cell-cell adhesion mediator activity. Involved in angiogenesis and establishment of endothelial barrier. Located in extracellular exosome. Implicated in aortic valve disease 3.

Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	MRB; AOVD3; ECSM4; ROBO4
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 848-1007 of human ROBO4 (NP_061928.4).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.