

# STRA13 Rabbit Polyclonal Antibody



CAB9892

---

## Product Information

### Size:

20uL, 50uL, 100uL, 200uL

### Observed MW:

14kDa

### Calculated MW:

6kDa/7kDa/8kDa

### Applications:

WB

### Reactivity:

Human

## Antibody Information

### Recommended dilutions:

WB 1:500 - 1:2000

### Source:

Rabbit

### Isotype:

IgG

### Purification:

Affinity purification

## Protein Background

### Immunogen information

#### Gene ID:

201254

#### Uniprot

A8MT69

#### Synonyms:

CENPX; CENP-X; D9; FAAP10; MHF2; STRA13

#### Immunogen:


Recombinant fusion protein containing a sequence corresponding to amino acids 1-63 of human STRA13 (NP\_659435.2).

#### Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---



Western blot analysis of extracts of H460 cells, using STRA13 antibody (CAB9892) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (CABS014) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit (CABM00020). Exposure time: 60s.