

Product Information

Product SKU:	CAB9894	Gene ID:	646486	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Rat

Additional Information

Observed MW:	16kDa	Conjugate:	Unconjugated
Calculated MW:	16kDa	Isotype:	IgG

Immunogen Information

Background: Predicted to enable lipid binding activity. Predicted to be located in cytosol.

Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	FABP12
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 1-140 of human FABP12 (NP_001098751.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.