

Product Information

Product SKU:	CAB9930	Gene ID:	971	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	40kDa	Conjugate:	Unconjugated
Calculated MW:	40kDa	Isotype:	IgG

Immunogen Information

**Background:** Predicted to enable signaling receptor binding activity. Predicted to be involved in cell adhesion. Is integral component of plasma membrane.

Recommended Dilution:	WB,1:500 - 1:1000
Synonyms:	LYB2; CD72b; CD72
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 120-359 of human CD72 (NP_001773.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.