

Product Information

Product SKU:	CAB9984	Gene ID:	80097	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human

Additional Information

Observed MW:	16kDa	Conjugate:	Unconjugated
Calculated MW:	16kDa	Isotype:	IgG

Immunogen Information

**Background:** Located in cytosol; microtubule cytoskeleton; and nucleoplasm. Part of gamma-tubulin large complex.

Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	FAM128B; MOZART2B; MZT2B
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 1-158 of human MZT2B (NP_079305.2).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.