

UBE2Q2 Rabbit Polyclonal Antibody

CAB9992



Product Information

Product SKU:	CAB9992	Gene ID:	92912	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse

Additional Information

Observed MW:	43kDa	Conjugate:	Unconjugated
Calculated MW:	43kDa	Isotype:	IgG

Immunogen Information

Background: Enables ubiquitin-protein transferase activity. Involved in protein K48-linked ubiquitination. Predicted to be located in cytosol. Predicted to be active in nucleus.

Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	UBE2Q2
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 1-150 of human UBE2Q2 (NP_775740.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.