

# Mouse Control IgG



CABC011

---

## Product Information

**Size:**

100 uL

**Observed MW:**

**Calculated MW:**

**Applications:**

IP

**Reactivity:**

## Protein Background

### Immunogen information

**Gene ID:**

Uniprot

**Synonyms:**

**Immunogen:**

None

## Antibody Information

**Recommended dilutions:**

IP 1:50 - 1:200

**Source:**

Mouse

**Isotype:**

IgG

**Purification:**


Protein A/G purification

**Storage:**

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---



Immunoprecipitation of over-expressed Myc-tagged protein in 293T cells incubated using Myc-tag antibody (CABE010). A mock served as negative control using mouse Control IgG (CABC011) and over-expressed 293T cell lysate served as positive control.