

Rat Control IgG



CABC034

Product Information

Size:

100 uL

Observed MW:

Calculated MW:

Applications:

Reactivity:

Antibody Information

Recommended dilutions:

Source:

Rat

Isotype:

IgG

Purification:

Protein A/G purification

Protein Background

Immunogen information

Gene ID:

Uniprot

Synonyms:

Immunogen:

None

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
