

CABP0093

Product Information

Product SKU:	CABP0093	Gene ID:	8290/8350	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat,Other (Wide Range Predicted)

Additional Information

Observed MW:	17kDa	Conjugate:	Unconjugated
Calculated MW:	15kDa	Isotype:	IgG

Immunogen Information

Background:	Histones are basic nuclear proteins that are responsible for the nucleosome structure of the chromosomal fiber in eukaryotes. This structure consists of approximately 146 bp of DNA wrapped around a nucleosome, an octamer composed of pairs of each of the four core histones (H2A, H2B, H3, and H4). The chromatin fiber is further compacted through the interaction of a linker histone, H1, with the DNA between the nucleosomes to form higher order chromatin structures. This gene is intronless and encodes a replication-dependent histone that is a member of the histone H3 family. Transcripts from this gene lack polyA tails; instead, they contain a palindromic termination element. This gene is found in the large histone gene cluster on chromosome 6p22-p21.3.
Recommended Dilution:	WB,1:1000 - 1:5000 IHC-P,1:50 - 1:200 IF/ICC,1:50 - 1:200
Synonyms:	H3/A; H3C2; H3C3; H3C4; H3C6; H3C7; H3C8; H3FA; H3C10; H3C11; H3C12; HIST1H3A; Phospho-Histone H3-T11
Purification Method:	Affinity purification
Immunogen:	A synthetic phosphorylated peptide around T11 of human HIST1H3B (NP_003520.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.