

CABP0105

Product Information

Product SKU:	CABP0105	Gene ID:	3725	Size:	20uL, 100uL
Clone No:	ARC0051	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	48kDa	Conjugate:	Unconjugated
Calculated MW:	36kDa	Isotype:	IgG

Immunogen Information

Background: This gene is the putative transforming gene of avian sarcoma virus 17. It encodes a protein which is highly similar to the viral protein, and which interacts directly with specific target DNA sequences to regulate gene expression. This gene is intronless and is mapped to 1p32-p31, a chromosomal region involved in both translocations and deletions in human malignancies.

Recommended Dilution: WB,1:500 - 1:2000 IF/ICC,1:50 - 1:200

Synonyms: AP1; p39; AP-1; cJUN; c-Jun; Phospho-c-Jun-S63

Purification Method: Affinity purification

Immunogen: A synthetic phosphorylated peptide around S63 of human c-Jun (P05412).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.