

# Phospho-JUNB-S79 Rabbit Polyclonal Antibody



CABP0376

---

## Product Information

### Size:

50uL, 100uL, 200uL

### Observed MW:

### Calculated MW:

35kDa

### Applications:

IHC IF

### Reactivity:

Human, Mouse, Rat

## Antibody Information

### Recommended dilutions:

IHC 1:50 - 1:100 IF 1:100 - 1:200

### Source:

Rabbit

### Isotype:

IgG

### Purification:

Affinity purification

## Protein Background

Transcription factor involved in regulating gene activity following the primary growth factor response. Binds to the DNA sequence 5'-TGA[CG]TCA-3'.

## Immunogen information

### Gene ID:

3726

### Uniprot

P17275

### Synonyms:

JUNB; AP-1

### Immunogen:


A phospho specific peptide corresponding to residues surrounding S79 of human JUNB

### Storage:


Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---



Immunohistochemistry of paraffin-embedded human breast carcinoma using Phospho-JUNB-S79 antibody (CABP0376).



Immunofluorescence analysis of methanol-fixed HeLa cells using Phospho-JUNB-S79 antibody (CABP0376).