

CABP0497

Product Information

Product SKU:	CABP0497	Gene ID:	5894	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human

Additional Information

Observed MW:	73kDa	Conjugate:	Unconjugated
Calculated MW:	73kDa	Isotype:	IgG

Immunogen Information

Background:	This gene is the cellular homolog of viral raf gene (v-raf). The encoded protein is a MAP kinase kinase kinase (MAP3K), which functions downstream of the Ras family of membrane associated GTPases to which it binds directly. Once activated, the cellular RAF1 protein can phosphorylate to activate the dual specificity protein kinases MEK1 and MEK2, which in turn phosphorylate to activate the serine/threonine specific protein kinases, ERK1 and ERK2. Activated ERKs are pleiotropic effectors of cell physiology and play an important role in the control of gene expression involved in the cell division cycle, apoptosis, cell differentiation and cell migration. Mutations in this gene are associated with Noonan syndrome 5 and LEOPARD syndrome 2.
Recommended Dilution:	WB, 1:500 - 1:2000 IP, 0.5µg-4µg antibody for 200µg-400µg extracts of whole cells
Synonyms:	NS5; CRAF; Raf-1; c-Raf; CMD1NN; Phospho-Raf1-S259
Purification Method:	Affinity purification
Immunogen:	A synthetic phosphorylated peptide around S259 of human Phospho-Raf1-S259 (NP_002871.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.