

CABP0534

Product Information

Product SKU:	CABP0534	Gene ID:	5587	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	102kDa	Conjugate:	Unconjugated
Calculated MW:	102kDa	Isotype:	IgG

Immunogen Information

Background:	The protein encoded by this gene is a serine/threonine protein kinase involved in many cellular processes, including Golgi body membrane integrity and transport, cell migration and differentiation, MAPK8/JNK1 and Ras pathway signaling, MAPK1/3 (ERK1/2) pathway signaling, cell survival, and regulation of cell shape and adhesion.
Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	PKD; PKCM; CHDED; PRKCM; PKC-MU; Phospho-PKC mu-S744/748
Purification Method:	Affinity purification
Immunogen:	A synthetic phosphorylated peptide around S744 & S748 of mouse PKC mu (NP_032884.2).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.