

Phospho-ULK1-S423 Rabbit Polyclonal Antibody



CABP0654

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Refer to figures

Calculated MW:

112kDa

Applications:

WB

Reactivity:

Human, Mouse, Rat

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

8408

Uniprot

O75385

Synonyms:

ULK1; ATG1; ATG1A; UNC51; Unc51.1; hATG1; ULK1

Immunogen:

A phospho specific peptide corresponding to residues surrounding S423 of human ULK1

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
