

CABP0684

Product Information

Product SKU:	CABP0684	Gene ID:	2309	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	97kDa	Conjugate:	Unconjugated
Calculated MW:	71kDa	Isotype:	IgG

Immunogen Information

Background: This gene belongs to the forkhead family of transcription factors which are characterized by a distinct forkhead domain. This gene likely functions as a trigger for apoptosis through expression of genes necessary for cell death. Translocation of this gene with the MLL gene is associated with secondary acute leukemia. Alternatively spliced transcript variants encoding the same protein have been observed.

Recommended Dilution: WB,1:500 - 1:2000 IHC-P,1:50 - 1:200

Synonyms: FOXO2; AF6q21; FKHRL1; FOXO3A; FKHRL1P2; Phospho-FOXO3A-S253

Purification Method: Affinity purification

Immunogen: A synthetic phosphorylated peptide around S253 of human FOXO3A (NP_001446.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.