

Phospho-eIF4B-S406 Rabbit Polyclonal Antibody



CABP0782

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Refer to Figures

Calculated MW:

64kDa/69kDa

Applications:

WB

Reactivity:

Human, Mouse, Rat

Protein Background

Required for the binding of mRNA to ribosomes. Functions in close association with EIF4-F and EIF4-A. Binds near the 5'-terminal cap of mRNA in presence of EIF-4F and ATP. Promotes the ATPase activity and the ATP-dependent RNA unwinding activity of both EIF4-A and EIF4-F.

Immunogen information

Gene ID:

1975

Uniprot

P23588

Synonyms:

EIF4B; EIF-4B; PRO1843

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Immunogen:

A synthetic phosphorylated peptide around S406 of human eIF4B (NP_001408.2).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
