

Phospho-H2AX-S134 Rabbit Polyclonal Antibody



CABP0888

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Refer to figures

Calculated MW:

14kDa

Applications:

WB

Reactivity:

Oryza sativa

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

4333939

Uniprot

Q8LLP5

Synonyms:

H2AX

Immunogen:

A phospho specific peptide corresponding to residues surrounding S134 of Oryza sativa H2AX.

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
