

Phospho-RALGAPA1-Thr735 Rabbit Polyclonal Antibody

CABP0908



Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

320kDa

Calculated MW:**Applications:**

WB IF

Reactivity:

Human, Mouse

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000 IF 1:50 - 1:200

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Catalytic subunit of the heterodimeric RalGAP1 complex which acts as a GTPase activator for the Ras-like small GTPases RALA and RALB (By similarity). May interact with the HLH region of TCF3/isoform E12.

Immunogen information

Gene ID:

56784

Uniprot

Q6GYP7

Synonyms:

2310003F20Rik; 4930400K19Rik; AI563624; GRIPE; Garnl1; Tulip1; mKIAA0884


Immunogen:

A phospho specific peptide corresponding to residues surrounding Thr735 of human RALGAPA1

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot - Phospho-RALGAPA1-Thr735 Rabbit pAb (CABP0908)