

Phospho-ERN1-Ser551 Rabbit Polyclonal Antibody

CABP0911



Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Refer to Figures

Calculated MW:

6kDa/109kDa

Applications:

WB

Reactivity:

Human, Rat

Protein Background

The protein encoded by this gene is the ER to nucleus signalling 1 protein, a human homologue of the yeast Ire1 gene product. This protein possesses intrinsic kinase activity and an endoribonuclease activity and it is important in altering gene expression as a response to endoplasmic reticulum-based stress signals.

Immunogen information

Gene ID:

2081

Uniprot

O75460

Synonyms:

ERN1; IRE1; IRE1P; IRE1a; hIRE1p

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Immunogen:

A phospho specific peptide corresponding to residues surrounding Ser551 of human ERN1.

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
