

Phospho-SMO-S615 Rabbit Polyclonal Antibody



CABP0937

Product Information

Size:

50uL, 100uL, 200uL

Observed MW:

Refer to figures

Calculated MW:

86kDa

Applications:

WB

Reactivity:

Human

Protein Background

The protein encoded by this gene is a G protein-coupled receptor that interacts with the patched protein, a receptor for hedgehog proteins. The encoded protein transduces signals to other proteins after activation by a hedgehog protein/patched protein complex.

Immunogen information

Gene ID:

6608

Uniprot

Q99835

Synonyms:

Gx; CRJS; SMOH; FZD11; SMO

Antibody Information

Recommended dilutions:

WB 1:500 - 1:1000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Immunogen:

A synthetic phosphorylated peptide around S615 of human SMO.

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
