

CABP0940

Product Information

Product SKU:	CABP0940	Gene ID:	6608	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human

Additional Information

Observed MW:	80kDa	Conjugate:	Unconjugated
Calculated MW:	86kDa	Isotype:	IgG

Immunogen Information

**Background:** The protein encoded by this gene is a G protein-coupled receptor that interacts with the patched protein, a receptor for hedgehog proteins. The encoded protein transduces signals to other proteins after activation by a hedgehog protein/patched protein complex.

Recommended Dilution:	WB, 1:500 - 1:1000
Synonyms:	Gx; CRJS; PHLS; SMOH; FZD11; Phospho-SMO-S611/615/616
Purification Method:	Affinity purification
Immunogen:	A synthetic phosphorylated peptide around S611 & S615 & S616 of human SMO (NP_005622.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.