

CABP0986

Product Information

Product SKU:	CABP0986	Gene ID:	7157	Size:	20uL, 100uL
Clone No:	ARC1530	Host Species:	Rabbit	Reactivity:	Human

Additional Information

Observed MW:	53kDa	Conjugate:	Unconjugated
Calculated MW:	44kDa	Isotype:	IgG

Immunogen Information

Background:	This gene encodes a tumor suppressor protein containing transcriptional activation, DNA binding, and oligomerization domains. The encoded protein responds to diverse cellular stresses to regulate expression of target genes, thereby inducing cell cycle arrest, apoptosis, senescence, DNA repair, or changes in metabolism. Mutations in this gene are associated with a variety of human cancers, including hereditary cancers such as Li-Fraumeni syndrome. Alternative splicing of this gene and the use of alternate promoters result in multiple transcript variants and isoforms. Additional isoforms have also been shown to result from the use of alternate translation initiation codons from identical transcript variants (PMIDs: 12032546, 20937277).
Recommended Dilution:	WB, 1:500 - 1:2000
Synonyms:	P53; BCC7; LFS1; BMFS5; TRP53; Phospho-p53-T55
Purification Method:	Affinity purification
Immunogen:	A synthetic phosphorylated peptide around T55 of human p53 (P04637).
Storage:	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 0.05% BSA, 50% glycerol, pH 7.3.