

Phospho-EGFR-Y1173 Rabbit Monoclonal Antibody

CABP0992



Product Information

Size:

20uL, 50uL, 100uL

Observed MW:

175KDa

Calculated MW:

175kDa

Applications:

WB

Reactivity:

Human

Protein Background

The protein encoded by this gene is a transmembrane glycoprotein that is a member of the protein kinase superfamily. This protein is a receptor for members of the epidermal growth factor family. EGFR is a cell surface protein that binds to epidermal growth factor. Binding of the protein to a ligand induces receptor dimerization and tyrosine autophosphorylation and leads to cell proliferation. Mutations in this gene are associated with lung cancer. [provided by RefSeq, Jun 2016]

Immunogen information

Gene ID:

1956

Uniprot

P00533

Synonyms:

ERBB; ERBB1; HER1; NISBD2; PIG61; mENA

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification


Immunogen:

A phospho specific peptide corresponding to residues surrounding Y1173 of human EGFR

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 0.05% BSA, 50% glycerol, pH7.3.

Product Images



Western blot - Phospho-EGFR-Y1173 Rabbit mAb (CABP0992)