

CABP0992

Product Information

Product SKU:	CABP0992	Gene ID:	1956	Size:	20uL, 100uL
Clone No:	ARC1536	Host Species:	Rabbit	Reactivity:	Human

Additional Information

Observed MW:	175kDa	Conjugate:	Unconjugated
Calculated MW:	134kDa	Isotype:	IgG

Immunogen Information

Background: The protein encoded by this gene is a transmembrane glycoprotein that is a member of the protein kinase superfamily. This protein is a receptor for members of the epidermal growth factor family. EGFR is a cell surface protein that binds to epidermal growth factor, thus inducing receptor dimerization and tyrosine autophosphorylation leading to cell proliferation. Mutations in this gene are associated with lung cancer. EGFR is a component of the cytokine storm which contributes to a severe form of Coronavirus Disease 2019 (COVID-19) resulting from infection with severe acute respiratory syndrome coronavirus-2 (SARS-CoV-2).

Recommended Dilution: WB, 1:500 - 1:1000

Synonyms: ERBB; ERBP; HER1; mENA; ERBB1; PIG61; NISBD2; Phospho-EGFR-Y1173

Purification Method: Affinity purification

Immunogen: A phospho specific peptide corresponding to residues surrounding Y1173 of human EGFR

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 0.05% BSA, 50% glycerol, pH 7.3.