

## CABP1003

---

**Product Information**

<b>Product SKU:</b>	CABP1003	<b>Gene ID:</b>	3725	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	ARC1548	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Human,Mouse,Rat

---

**Additional Information**

<b>Observed MW:</b>	43kDa	<b>Conjugate:</b>	Unconjugated
<b>Calculated MW:</b>	36kDa	<b>Isotype:</b>	IgG

---

**Immunogen Information**

**Background:** This gene is the putative transforming gene of avian sarcoma virus 17. It encodes a protein which is highly similar to the viral protein, and which interacts directly with specific target DNA sequences to regulate gene expression. This gene is intronless and is mapped to 1p32-p31, a chromosomal region involved in both translocations and deletions in human malignancies.

**Recommended Dilution:** WB,1:500 - 1:2000 IHC-P,1:50 - 1:200

**Synonyms:** AP1; p39; AP-1; cJUN; c-Jun; Phospho-c-Jun-T91

**Purification Method:** Affinity purification

**Immunogen:** A synthetic phosphorylated peptide around T91 of human c-Jun (P05412).

**Storage:** Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.