

Phospho-c-Jun-T91 Rabbit Monoclonal Antibody



CABP1003

Product Information

Size:

20uL, 50uL, 100uL

Observed MW:

43KDa

Calculated MW:

43kDa

Applications:

WB IHC

Reactivity:

Human, Mouse, Rat

Protein Background

This gene is the putative transforming gene of avian sarcoma virus 17. It encodes a protein which is highly similar to the viral protein, and which interacts directly with specific target DNA sequences to regulate gene expression. This gene is intronless and is mapped to 1p32-p31, a chromosomal region involved in both translocations and deletions in human malignancies. [provided by RefSeq, Jul 2008]

Immunogen information

Gene ID:

3725

Uniprot

P05412

Synonyms:

AP-1; AP1; c-Jun

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000 IHC 1:50
- 1:200

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification


Immunogen:

A phospho specific peptide corresponding to residues surrounding T91 of human c-Jun

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 0.05% BSA, 50% glycerol, pH7.3.

Product Images



Western blot - Phospho-c-Jun-T91 Rabbit mAb (CABP1003)