

CABP1020

Product Information

Product SKU:	CABP1020	Gene ID:	5604	Size:	20uL, 100uL
Clone No:	ARC1565	Host Species:	Rabbit	Reactivity:	Human,Rat

Additional Information

Observed MW:	45kDa	Conjugate:	Unconjugated
Calculated MW:	43kDa	Isotype:	IgG

Immunogen Information

Background: The protein encoded by this gene is a member of the dual specificity protein kinase family, which acts as a mitogen-activated protein (MAP) kinase kinase. MAP kinases, also known as extracellular signal-regulated kinases (ERKs), act as an integration point for multiple biochemical signals. This protein kinase lies upstream of MAP kinases and stimulates the enzymatic activity of MAP kinases upon wide variety of extra- and intracellular signals. As an essential component of MAP kinase signal transduction pathway, this kinase is involved in many cellular processes such as proliferation, differentiation, transcription regulation and development.

Recommended Dilution: WB,1:500 - 1:1000

Synonyms: MEL; Phospho-MEK1-S298

Purification Method: Affinity purification

Immunogen: A synthetic phosphorylated peptide around S298 of human Phospho-MEK1-S298 (Q02750).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.