

Phospho-Lyn-Y396 Rabbit Monoclonal Antibody



CABP1050

Product Information

Size:

20uL, 50uL, 100uL

Observed MW:

56KDa

Calculated MW:

56kDa

Applications:

WB

Reactivity:

Human, Mouse, Rat

Protein Background

This gene encodes a tyrosine protein kinase, which maybe involved in the regulation of mast cell degranulation, and erythroid differentiation. Alternatively spliced transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2011]

Immunogen information

Gene ID:

4067

Uniprot

P07948

Synonyms:

JTK8; p53Lyn; p56Lyn

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification


Immunogen:

A phospho specific peptide corresponding to residues surrounding Y396 of human Lyn

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 0.05% BSA, 50% glycerol, pH7.3.

Product Images



Western blot - Phospho-Lyn-Y396 Rabbit mAb (CABP1050)