

## CABP1050

---

**Product Information**

<b>Product SKU:</b>	CABP1050	<b>Gene ID:</b>	4067	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	ARC1597	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Human,Mouse,Rat

---

**Additional Information**

<b>Observed MW:</b>	56kDa	<b>Conjugate:</b>	Unconjugated
<b>Calculated MW:</b>	56kDa/59kDa	<b>Isotype:</b>	IgG

---

**Immunogen Information**

**Background:** This gene encodes a tyrosine protein kinase, which maybe involved in the regulation of mast cell degranulation, and erythroid differentiation. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.

**Recommended Dilution:** WB,1:500 - 1:2000

**Synonyms:** JTK8; p53Lyn; p56Lyn; Phospho-Lyn-Y396

**Purification Method:** Affinity purification

**Immunogen:** A synthetic phosphorylated peptide around Y396 of human Lyn (P07948).

**Storage:** Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.