

Phospho-Delta-Catenin/p120 Catenin-S320 Rabbit Polyclonal Antibody CABP1079



Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Refer to figures

Calculated MW:

65-72kDa/90-108kDa

Applications:

WB

Reactivity:

Human

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

This gene encodes a member of the Armadillo protein family, which function in adhesion between cells and signal transduction. Multiple translation initiation codons and alternative splicing result in many different isoforms being translated. Not all of the full-length nature of the described transcript variants have been determined. Read-through transcription also exists between this gene and the neighboring upstream thioredoxin-related transmembrane protein 2 (TMX2) gene. [provided by RefSeq, Dec 2010]

Immunogen information

Gene ID:

1500

Uniprot

O60716

Synonyms:

CAS; CTNND; P120CAS; P120CTN; p120; p120(CAS); p120(CTN); CTNND1

Immunogen:

A phospho synthetic peptide corresponding to residues surrounding S320 of human Delta-Catenin/p120 Catenin.

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
