

Phospho-STK11-S428 Rabbit Polyclonal Antibody



CABP1108

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

65kDa

Calculated MW:

45kDa/48kDa

Applications:

WB

Reactivity:

Human, Mouse, Rat

Protein Background

This gene, which encodes a member of the serine/threonine kinase family, regulates cell polarity and functions as a tumor suppressor. Mutations in this gene have been associated with Peutz-Jeghers syndrome, an autosomal dominant disorder characterized by the growth of polyps in the gastrointestinal tract, pigmented macules on the skin and mouth, and other neoplasms. Alternate transcriptional splice variants of this gene have been observed but have not been thoroughly characterized. [provided by RefSeq, Jul 2008]

Immunogen information

Gene ID:

6794

Uniprot

Q15831

Synonyms:

LKB1; PJS; hLKB1; STK11

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification


Immunogen:

A phospho specific peptide corresponding to residues surrounding S428 of human STK11.

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot - Phospho-STK11-S428 Rabbit pAb (CABP1108)