

Phospho-DAPP1-Y139 Rabbit Polyclonal Antibody

CABP1128



Product Information

Size:

50uL, 100uL, 200uL

Observed MW:

Refer to figures

Calculated MW:

32kDa

Applications:

WB

Reactivity:

Human, Mouse

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

May act as a B-cell-associated adapter that regulates B-cell antigen receptor (BCR)-signaling downstream of PI3K.

Immunogen information

Gene ID:

27071

Uniprot

Q9UN19

Synonyms:

DAPP1; BAM32; dual adaptor of phosphotyrosine and 3-phosphoinositides 1

Immunogen:

A phospho specific peptide corresponding to residues surrounding Y139 of human DAPP1.

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
