

# Phospho-FOXO3-S253 Rabbit Monoclonal Antibody

## CABP1131



### Product Information

**Size:**

20uL, 50uL, 100uL

**Observed MW:****Calculated MW:**

97kDa

**Applications:**

WB

**Reactivity:**

Human, Mouse, Rat

### Protein Background

This gene belongs to the forkhead family of transcription factors which are characterized by a distinct forkhead domain. This gene likely functions as a trigger for apoptosis through expression of genes necessary for cell death. Translocation of this gene with the MLL gene is associated with secondary acute leukemia. Alternatively spliced transcript variants encoding the same protein have been observed. [provided by RefSeq, Jul 2008]

### Immunogen information

**Gene ID:**

2309

**Uniprot**

O43524

**Synonyms:**

AF6q21; FKHL1; FKHL1P2; FOXO2; FOXO3A

### Antibody Information

**Recommended dilutions:**

WB 1:500 - 1:2000

**Source:**

Rabbit

**Isotype:**

IgG

**Purification:**

Affinity purification

**Immunogen:**

A phospho specific peptide corresponding to residues surrounding S253 of human FOXO3

**Storage:**

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 0.05% BSA, 50% glycerol, pH7.3.

## Product Images

---