

Phospho-Paxillin-Y118 Rabbit Monoclonal Antibody

CABP1156



Product Information

Size:

20uL, 50uL, 100uL

Observed MW:

75kDa

Calculated MW:

68kDa

Applications:

WB

Reactivity:

Human, Mouse, Rat

Protein Background

This gene encodes a cytoskeletal protein involved in actin-membrane attachment at sites of cell adhesion to the extracellular matrix (focal adhesion). Alternatively spliced transcript variants encoding different isoforms have been described for this gene. These isoforms exhibit different expression pattern, and have different biochemical, as well as physiological properties (PMID:9054445). [provided by RefSeq, Aug 2011]

Immunogen information

Gene ID:

5829

Uniprot

P49023

Synonyms:

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Immunogen:

A phospho specific peptide corresponding to residues surrounding Y118 of human Paxillin

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 0.05% BSA, 50% glycerol, pH7.3.

Product Images
