

HRP Goat Anti-Swine IgG (H+L)



CABS004

Product Information

Size:

100 uL

Observed MW:

Calculated MW:

Applications:

WB IHC

Reactivity:

Protein Background

Immunogen information

Gene ID:

Uniprot

Synonyms:

Immunogen:

Swine IgG

Antibody Information

Recommended dilutions:

WB 1:2000 - 1:10000 IHC
1:100 - 1:1000 ICC 1:100 -
1:1000

Source:

Goat

Isotype:

Horseradish peroxidase
conjugated IgG

Purification:

Affinity purification

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 50%
glycerol, pH7.3.

Product Images
