

CABS014

Product Information

Product SKU:	CABS014	Gene ID:	-	Size:	20uL, 100uL
Clone No:	-	Host Species:	Goat	Reactivity:	-

Additional Information

Observed MW:	-	Conjugate:	HRP
Calculated MW:	-	Isotype:	Horseradish peroxidase conjugated IgG

Immunogen Information

Background:	Secondary antibodies are affinity-purified antibodies which will work with target-specific primary antibody in the detection, sorting or purification of its specified target. Secondary antibodies offer increased versatility enabling users to use many detection systems (e.g. HRP, AP, fluorescence). They can also provide greater sensitivity through signal amplification as multiple secondary antibodies. Most commonly, secondary antibodies are generated by immunizing the host animal (different from host species of primary antibody) with a pooled population of normal immunoglobulins from the host species of primary antibody and can be further purified and modified (i.e. antibody fragmentation, label conjugation, etc.) to ensure well-characterized specificity to corresponding normal immunoglobulins.
Recommended Dilution:	ELISA, 1:5000 - 1:10000 WB, 1:2000 - 1:10000 IHC-P, 1:50 - 1:200 DB, 1:2000 - 1:10000
Synonyms:	-
Purification Method:	Affinity purification
Immunogen:	Rabbit IgG
Storage:	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA, 50% glycerol, pH 7.3.