

TRITC Goat Anti-Mouse IgG (H+L)



CABS026

Product Information

Size:

100 uL

Observed MW:

Calculated MW:

Applications:

IF FC

Reactivity:

Protein Background

Immunogen information

Gene ID:

Uniprot

Synonyms:

Immunogen:

Mouse IgG

Antibody Information

Recommended dilutions:

IF 1:50 - 1:200 FC 1:50 -
1:200

Source:

Goat

Isotype:

TRITC conjugated IgG

Purification:

Affinity purification

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: 50% glycerol,
pH7.3.

Product Images
