

HRP Rabbit Anti-Goat IgG (H+L)



CABS029

Product Information

Size:

100 uL

Observed MW:

Calculated MW:

Applications:

WB IHC

Reactivity:

Protein Background

Immunogen information

Gene ID:

Uniprot

Synonyms:

Immunogen:

Goat IgG

Antibody Information

Recommended dilutions:

WB 1:5000 - 1:100000 IHC

1:500 - 1:5000 ICC 1:500 -

1:5000

Source:

Rabbit

Isotype:

Horseradish peroxidase

conjugated IgG

Purification:

Affinity purification

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: 50% glycerol, pH7.3.

Product Images
