

CABS040

Product Information

Product SKU:	CABS040	Gene ID:	-	Size:	20uL, 100uL
Clone No:	-	Host Species:	Goat	Reactivity:	-

Additional Information

Observed MW:	-	Conjugate:	Rhodamine. Ex:550nm. Em:570nm.
Calculated MW:	-	Isotype:	TRITC conjugated IgG

Immunogen Information

Background:	Secondary antibodies are affinity-purified antibodies which will work with target-specific primary antibody in the detection, sorting or purification of its specified target. Secondary antibodies offer increased versatility enabling users to use many detection systems (e.g. HRP, AP, fluorescence). They can also provide greater sensitivity through signal amplification as multiple secondary antibodies. Most commonly, secondary antibodies are generated by immunizing the host animal (different from host species of primary antibody) with a pooled population of normal immunoglobulins from the host species of primary antibody and can be further purified and modified (i.e. antibody fragmentation, label conjugation, etc.) to ensure well-characterized specificity to corresponding normal immunoglobulins.
Recommended Dilution:	IHC-P, 1:50 - 1:200 IF/ICC, 1:50 - 1:200 FC, 1:50 - 1:200
Synonyms:	-
Purification Method:	Affinity purification
Immunogen:	Rabbit IgG
Storage:	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.025% Sodium Azide, 0.75% BSA, 50% glycerol, pH 7.3.