

CABS058

Product Information

Product SKU:	CABS058	Gene ID:	-	Size:	20uL, 100uL
Clone No:	-	Host Species:	Goat	Reactivity:	Wide range of species

Additional Information

Observed MW:	-	Conjugate:	Alexa Fluor 555. Ex:555nm. Em:565nm.
Calculated MW:	-	Isotype:	IgG

Immunogen Information

Background:	Secondary antibodies are affinity-purified antibodies which will work with target-specific primary antibody in the detection, sorting or purification of its specified target. Secondary antibodies offer increased versatility enabling users to use many detection systems (e.g. HRP, AP, fluorescence). They can also provide greater sensitivity through signal amplification as multiple secondary antibodies . Most commonly, secondary antibodies are generated by immunizing the host animal (different from host species of primary antibody) with a pooled population of normal immunoglobulins from the host species of primary antibody and can be further purified and modified (i.e. antibody fragmentation, label conjugation, etc.) to ensure well-characterized specificity to corresponding normal immunoglobulins.
Recommended Dilution:	IF/ICC,1:100 - 1:500 FC,1:500 - 1:2000
Synonyms:	-
Purification Method:	Affinity purification
Immunogen:	Rabbit IgG
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.