

**Alexa Fluor 647-conjugated Goat Anti-Mouse IgG  
(H+L)  
CABS059**



---

**Product Information**

**Size:**

100 uL

**Observed MW:**

**Calculated MW:**

**Applications:**

IF FC

**Reactivity:**

**Protein Background**

**Immunogen information**

**Gene ID:**

Uniprot

**Synonyms:**

**Immunogen:**

Mouse IgG

**Antibody Information**

**Recommended dilutions:**

IF 1:100 - 1:500 FC 1:500 -  
1:2000

**Source:**

Goat

**Isotype:**

IgG

**Purification:**

Affinity purification

**Storage:**

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02%  
sodium azide, 50% glycerol, pH7.3.

## Product Images

---