

Alexa Fluor 647-conjugated AffiniPure Goat Anti-Mouse IgG (H+L)

CABS078



Product Information

Size:

100 uL

Observed MW:

Calculated MW:

Applications:

IF, FC

Reactivity:

Protein Background

Immunogen information

Gene ID:

Uniprot

Synonyms:

Immunogen:

Mouse IgG

Antibody Information

Recommended dilutions:

IF 1:100 - 1:800

FC 1:100 - 1:800

Source:

Goat

Isotype:

Alexa Fluor 647
conjugated IgG

Purification:

Affinity purification

Storage:

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.