

## CB6 Biosimilar (Anti-SARS-CoV-2) Antibody

### Product Information

**Product SKU:**

HDBS0011

**Size:**

100 µg

**Isotype:**

IgG1

**Host:**

Homo sapiens

**Category:**

Biosimilar Antibody

**Uniprot:**

P0DTC2

**Reactivity:**

SARS-CoV-2

**Synonyms:**

S protein RBD, Spike glycoprotein Receptor-binding domain, S glycoprotein RBD, Spike protein RBD, COVID-19

**Common Name:**

CB6

**Target:**

SARS

**Formulation:**

Powder

**Buffer:**

Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions.

**Guarantee:**

12 months from date of dispatch

**Storage and Stability:**

Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized antibodies are shipped at ambient temperature.

**Applications**

ELISA

**Recommended Dilutions:**

ELISA 1:5000-10000

**Research Use Only:**

Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.