

## Sacituzumab Biosimilar (Anti-Trop2) Antibody

### Product Information

**Product SKU:**

HDBS0023

**Size:**

100 µg

**Isotype:**

IgG1

**Host:**

Humanized

**Category:**

Biosimilar Antibody

**Uniprot:**

P09758

**Reactivity:**

Human

**Synonyms:**

TACSTD2,GA733-1,M1S1,TROP2

**Common Name:**

IMMU-132, hRS7-SN-38, hRS7-SN-38-ADC, hRS7-[CL-SN-38], hrS7-SN-38

**Target:**

Trop2

**Formulation:**

Powder

**Buffer:**

Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions.

**Guarantee:**

12 months from date of dispatch

**Storage and Stability:**

Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized antibodies are shipped at ambient temperature.

**Applications**

ELISA, Flow Cyt

**Recommended Dilutions:**

ELISA 1:5000-10000; Flow Cyt 1:100

**Research Use Only:**

Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

