

# Tigatuzumab Biosimilar (Anti-TNFRSF10B) Antibody

## Product Information

**Product SKU:**

HDBS0032

**Size:**

100 µg

**Isotype:**

IgG1

**Host:**

Humanized

**Category:**

Biosimilar Antibody

**Uniprot:**

O14763

**Reactivity:**

Human

**Synonyms:**

TNFRSF10B, TRAILR2, TRAIL-R2, CD262, DR5, KILLER, TRICK2, ZTNFR9, TRICKB

**Common Name:**

GEN-1029, HexaBody-DR5/DR5, Hx-DR5-01/05

**Target:**

TNFRSF10B

**Formulation:**

Powder

**Buffer:**



Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions.

**Guarantee:**

12 months from date of dispatch

**Storage and Stability:**

Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized antibodies are shipped at ambient temperature.

**Applications**

ELISA

**Recommended Dilutions:**

ELISA 1:5000-10000

**Research Use Only:**

Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

