

## Labetuzumab Biosimilar (Anti-CEACAM5) Antibody

### Product Information

**Product SKU:**

HDBS0035

**Size:**

100 µg

**Isotype:**

IgG1

**Host:**

Humanized

**Category:**

Biosimilar Antibody

**Uniprot:**

P06731

**Reactivity:**

Human

**Synonyms:**

CEACAM-5, CD66e, CEA, Meconium antigen 100

**Common Name:**

IMMU-130, hMN-14-SN-38, hMN-14-SN-38 ADC, hMN14-CL-SN-38, hMN14-SN-38

**Target:**

CEACAM5

**Formulation:**

Powder

**Buffer:**

Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions.



**Guarantee:**

12 months from date of dispatch

**Storage and Stability:**

Store at  $-20^{\circ}\text{C}$  to  $-80^{\circ}\text{C}$  for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at  $-80^{\circ}\text{C}$  (Avoid repeated freezing and thawing). Lyophilized antibodies are shipped at ambient temperature.

**Applications**

ELISA

**Recommended Dilutions:**

ELISA 1:5000-10000

**Research Use Only:**

Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

