

Ruplizumab Biosimilar (Anti-CD40L) Antibody

Product Information

Product SKU:

HDBS0042

Size:

100 µg

Isotype:

IgG1

Host:

Humanized

Category:

Biosimilar Antibody

Uniprot:

P29965

Reactivity:

Human

Synonyms:

CD40LG,CD154,CD40L,HIGM1,IGM,IMD3,T-BAM,TNFSF5,TRAP,gp39

Common Name:

BG-9588, hu5c8

Target:

CD40 Ligand

Formulation:

Powder

Buffer:

Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions.



Guarantee:

12 months from date of dispatch

Storage and Stability:

Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized antibodies are shipped at ambient temperature.

Applications

ELISA, Flow Cyt

Recommended Dilutions:

ELISA 1:5000-10000; Flow Cyt 1:100

Research Use Only:

Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

